中国人民大学化学与生命资源学院
SCHOOL OF CHEMISTRY AND LIFE RESOURCES, RENMIN UNIVERSITY OF CHINA

理化分析测试中心
INSTRUMENTAL ANALYSIS CENTER (IAC)

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1 Choose Cleaning Method(s)

Cell Cleaning Method	Syringe Cleaning Method
<input type="radio"/> Rinse Rinse with water.	<input checked="" type="radio"/> Rinse Rinse with water, then dry using methanol.
<input checked="" type="radio"/> Wash Wash with detergent, then rinse with water.	<input type="radio"/> Wash Wash with detergent, rinse with water, then dry using methanol.
<input type="radio"/> Soak Soak in detergent for 30 minutes at 60 °C, then rinse with water.	<input type="radio"/> None
<input type="radio"/> None	

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Malvern MicroCal PEAQ-ITC Control Software

Run Experiment Maintenance Design Experiment

Run Start Experiment Load Run Clean X

Experiment Information

[Syr] (M)	1.00e-3
[Cell] (M)	100e-6

Comment

Instrument Settings

Temperature (°C) 25.0

This experimental method will make a single 0.4 µL injection, followed by 18, 2 µL injections.


The following settings will be used:

Reference Power (µcal/s)	5.00
Feedback	High
Stir Speed (rpm)	750
Initial Delay (s)	60
Injection Spacing (s)	150
Injection Duration (s)	4.0

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


1 Load Cell

Move the pipette out of the way (i.e to the Clean Location).

Fill the loading syringe with 300 μ l of sample.

Slowly insert the loading syringe into the sample cell port, gently touch the cell bottom, and move up approximately 1 mm.

Slowly dispense approximately 150 μ l of sample.

Introduce a small volume, quickly, to dislodge any bubbles. Repeat this several times.

Slowly dispense the remaining sample while being careful not to introduce air bubbles.

Remove any excess sample in the cell port overflow cup using the loading syringe.

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2 Attach Fill Port Adapter

If the pipette is in the Clean Location, you must press the clamp's release lever.

Move the pipette to the Rest Location.

Align the hole in the pipette's housing to the hole in the pipette's rotating assembly.

Insert the fill port adapter. A soft click should be felt.

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3 Move Pipette to Load Location

Load approximately 60 μ l of titrant in one of the supplied microcentrifuge tubes.

Ensure the microcentrifuge tube has its lid properly situated in the keyed Load Location.

Move the pipette to the Load Location.

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4 Detach Fill Port Adapter

Move the pipette to the Rest Location.

Detach the fill port adapter from the pipette and return it to its Storage Location.

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